

## MAIN APPLICATIONS

Emergency Dept.

Home-Care

Patient Transfer

Room / Ward

Tracheostomy

Veterinary

Internal Rechargeable Battery



NEW ASKIR 230/12V BR is an electric suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. Suitable for emergency applications and for post-operative therapy, both in professional environment and home-care.

Three different options for operation: AC/DC adapter, rechargeable battery and 12V make it a versatile suction unit approved for NON-STOP operation without overheating. Provided with visual and acoustic alarm indicating low battery level.

The newly-designed ergonomic handle ensures a comfortable and secure grip.

## ACCESSORIES INCLUDED

- Liquid Collection Jar with overflow valve system (different options, see below)
- Silicone Tubes  $\varnothing$  6x10mm (autoclavable) - Patient Tube length 140cm
- Conical Connector  $\varnothing$  8-9-10mm
- Antibacterial & Hydrophobic Filter (single-patient) + CH20 Canula
- AC/DC Universal Adapter + Power Cord with Schuko plug
- 12V Car Adapter

## OPTIONS

- Stand on 5 Castors  $\varnothing$  50mm with Brakes (antistatic upon request) REF 27730
- Carrying Bag (for 1000ml jar) SP 0207/02
- Carrying Bag (for 2000ml jar) SP 0207
- Wall Mount 27736

## TECHNICAL FEATURES

Motor	Oilless and maintenance-free piston pump
Power Feeding	14V $\rightleftharpoons$ 4A by AC/DC Universal Adapter 100-240V / 50-60Hz - 100VA (model UE60-140429SPA1)
	12V $\rightleftharpoons$ 4A Internal rechargeable Pb battery - Autonomy 80 min. - Recharging Time 4 hours
	12V $\rightleftharpoons$ 1,9A Car adapter
ISO 10079-1 Classification	HIGH VACUUM / LOW FLOW
IP Code	IP21
Max Vacuum (adjustable)	-0.75 bar -75 kPa -563 mmHg
Max free air flow rate	16 l/min
Noise Level	63,0 dBA
Duty cycle	Non-stop operation (by mains only)
Weight	kg 4.55 (unit packed with all accessories)
Size	cm 35 x 18 x 21
Years of Warranty	2
Shipping carton	3
Place of Manufacturing	Italy

Italian  
Medical  
Touch  
MADE IN ITALY

## AVAILABLE CONFIGURATIONS

RE 310211	NEW ASKIR 230/12V BR with 1000 ml autoclavable jar in polycarbonate Makrolon 2858
RE 310211/01	NEW ASKIR 230/12V BR with 2000 ml autoclavable jar in polycarbonate Makrolon 2858
RE 310211/03	NEW ASKIR 230/12V BR with 1000 ml Flovac® disposable liner and reusable container
RE 310211/04	NEW ASKIR 230/12V BR with 2000 ml Flovac® disposable liner and reusable container
	Jars in APEC® autoclavable up to 143°C available upon request
	Disposable jar 1000 ml available upon request

CE 0123