



FAZZINI

Your Health is Our Aim

Fazzini srl
S.S Padana Superiore, 317-20055
Vimodrone-Milano-ITALY
Tel. 022651521 Fax. 0226515221
info@fazzini.it - www.fazzini.it

MECHANICAL SCALE

S7200EU



DIRECTIVE 2014/31 UE
DIRECTIVE EN45501 8.2

DIMENSIONS:

Length mm 530

Depth mm 280

Height mm 1510

Weight Kg 15

Base size mm 375x275

TECHNICAL FEATURES:

Mechanical column scale with eye-level beam. Large anti-slip feet surface.
4 adjustable feet and 2 castors for easy movement.

Base: with castors

Max. Capacity: 200 kg

Min. Capacity: 2 kg

Graduation: 100 gr

Tare: 2,2 kg