

panalex
THE O.R. COMPANY



SHANGHAI PAX MEDICAL INSTRUMENT CO., LTD.

Add: No. 3 Building No.438 Fulian Two Road.,Baoshan District Shanghai, China
P.C.:201906

Tel: +86-13601688868
Fax: +86-21-52820763
info@pax-medical.com
www.pax-medical.com

DS-IIM
Electric Operating Table



DS-II(M)

Electric Operating Table



DS-II(M) Electric Operating Table

- DS-II(M) electric operating table is designed for wide variety of surgical operations, including general operation, obstetrics, urology, ENT, endoscopic, neurosurgery, orthopedics, etc.
- There are standard table top and sliding table top available
- Wide sliding range of table top: 350mm, which can fully meet the requirement of X-ray C-arm
- Table frame is made of high strength cast aluminum alloy and stainless steel
- Memory foam mattress—comfortable for patients and prevent pressure ulcers
- 24V low voltage and automatic electric lock guarantee the operation safety
- Battery charge ensures the table can be used when power failure
- Automatic 0-position with one key for emergency case
- Side rail with European or USA standard that can truly meet the international market needs

Electric Operating Table

Loading weight	250kg
Sizes	55 cm x 210 cm
Adjustable height	65 - 110 cm (measured on the top of mattress)
Horizontal movement Stroke:	≥350mm (for C-ARM or X-RAY image taken)
Head section	Up>30° Down>55°, detachable
Left / right angle	20°
Reversed Trandelenburg	30°
Trandelenburg	25°
Kidney bridge elevation:	≥110mm
Leg section	down≥90°/ open≥180°, detachable



Independent kidney bridge design



Wide range of accessories for various operations



The two simultaneous control systems, hand control and control panel in column, ensure the table can work reliably in case one system down



PANALEX LED operating light